
ANEROID SPHYGMOMANOMETERS

seca b40/b41
seca b30/b31
seca b10/b11/b12

USER MANUAL
ANEROID SPHYGMOMANOMETERS

English	4–11
Français	12–17
Deutsch	18–23
Español	24–29
Suomi	30–35
Polski	36–41

EXPLANATION OF SYMBOLS | LÈGENDE | ERKLÄRUNG DER SYMBOLE | EXPLICACIÒN DE SÌMBOLOS | SYMBOLIEN SELITYKSET | OBJAŚNIENIE SYMBOLI:

Follow Instructions for Use.
Observer les consignes
du manuel de l'utilisateur.
Gebrauchsanweisung beachten.
Tener en cuenta las
instrucciones de uso.
Noudata käyttöohjetta.
Przestrzegać instrukcji obsługi.

Manufacturer
Fabricant
Hersteller
Fabricante
Valmistaja
Producent

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This product complies with the Council Directive 93/42/EEC concerning medical devices.

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„MR Conditional“ denotes suitability for use in a magnetic resonance (MR) environment in accordance with requirements ASTM F2052 – 06e1 (Standard Test Method For Measurement Of Magnetically Induced Displacement Force On Medical Devices In Magnetic Resonance Environment).

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TABLE OF CONTENTS

1. Description of device	6
2. Safety precautions	7
3. Operation	8
4. Hygiene treatment	9
5. Servicing	10
6. Disposal	10
7. Warranty	11
8. Declaration of conformity	11

1. DEVICE DESCRIPTION

1.1 Congratulations!

With this blood pressure measuring device, you have just purchased a precise and simultaneously robust device.

seca has used its technical expertise in the service of health care for more than 170 years and, as a market leader, has always set standards worldwide with its innovative developments.

1.2 Intended use

seca aneroid blood pressure measuring devices are used to measure blood pressure by auscultation on healthy skin on the upper arm or thigh.

1.3 Measuring range

The read-off scale has a measuring range from 0 – 300 mmHg.
Measuring precision across the whole measuring range is +/- 3 mmHg.

1.4 Contraindications

No contraindications are known to date.

2. SAFETY PRECAUTIONS

2.1 Safety precautions in these Instructions for Use

CAUTION!

Identifies a hazardous situation. If you fail to take note of this information, minor to moderate injury may result.

NOTICE!

Indicates that the product may have been operated incorrectly. If you fail to take note of this information, the device may be damaged or the measured results may be incorrect.

NOTE

Additional information on how to use this device.

2.2 Basic safety precautions on handling the device

- Please follow the precautions in these Instructions for Use.
- Keep the instructions for use in a safe place.

CAUTION! Avoiding injuries

- Young children have lower blood pressure than adults, do not pump up excessively.
- Measuring time should not exceed 2 minutes.
- Either use the cuff with protective hygiene strips or clean it before using it on the next patient.

NOTICE!

- Check that the product is working perfectly before using it each time.
- Do not use the device if it is visibly damaged.
- Under no circumstances pump up the pressure gage beyond 300 mmHg.
- Except for the seca b11 and seca b12, the products may not be taken into strong magnetic fields, for example MRI, and used there.
- Only use genuine seca spare parts.
- Avoid subjecting the device to shocks.
- Follow the care instructions.
- Avoid damaging the rubber parts with sharp or pointed objects.
- Do not expose the cuff and the device to any direct UV radiation, this applies to all plastic parts in particular.

NOTE

- Only qualified persons may interpret the measured results.

2.3 Important for domestic users

Consult your doctor regularly. He will let you know your individual normal blood pressure values and the value above which your blood pressure level has to be classified dangerous.

2.4 Ambient conditions

Ideal storage/use temperatures are between 0 °C and +40 °C with a relative humidity of below 85%.

3. OPERATION

3.1 Selecting cuff size

First select a suitably-sized cuff for the circumference of the extremity. The table below provides information about the relevant sphere of application.

Size 1	Infant	10.0 – 15.0 cm circumference	Free from PVC, DEHP and latex
Size 2	Young child	14.0 – 21.5 cm circumference	Free from PVC, DEHP and latex
Size 3	Adult/small	20.5 – 28.0 cm circumference	Free from PVC, DEHP and latex
Size 4	Adult	27.0 – 35.0 cm circumference	Free from PVC, DEHP and latex
Size 5	Adult/large	34.0 – 43.0 cm circumference	Free from PVC, DEHP and latex
Size 6	Thigh	42.0 – 54.0 cm circumference	Free from PVC, DEHP and latex

3.2 Performing a measurement

The cuff is applied so that the lower edge of the cuff is approx. 2 to 3 cm above the elbow/approx. 5 cm above the knee joint. Note the marking for artery position as you do this. Wind the cuff around the extremity neither too loosely nor too tightly, and ensure that the edge of the marking is within the permitted range. Close the valve on the air drain screw. The pressure required (approx. 30 mmHg above expected systolic pressure) is generated in the cuff by pumping the inflation bulb just a few times. The valve is opened slowly and blow-off speed metered to 2 to 3 mmHg per second. Using a stethoscope, it is now possible to auscultate Korotkoff sounds. Once measurement is complete, open the valve to vent the cuff completely.

NOTE

- Modification for left and right-handed use, seca b10, b11 and b12

1.

2.

- MR suitability of the seca b11 and b12

The seca b11 and b12 pressure gages were tested for suitability for use in a magnetic resonance (MR) environment in combination with the cuffs and their Velcro closures. Both the requirements on functionality and on electromagnetic attractive force were met with a range using 3 Tesla. This is indicated by the marking "MR Conditional" on the scale.

4. HYGIENE TREATMENT

4.1 Cleaning and caring for the pressure gage

NOTICE! Damage to device

Unsuitable cleaning agents and disinfectants may damage the sensitive surfaces of the device.

- Do not use aggressive or abrasive cleaning agents.
- Do not use white spirit or petroleum spirit.
- Only use disinfectants suitable for sensitive surfaces. Suitable disinfectants are available from specialist dealers.
- For cleaning, we recommend a damp cotton or microfiber cloth wetted with disinfectant if appropriate.
- Prevent liquids from penetrating the housing.

4.2 Cleaning and caring for the cuffs

Care information for seca blood pressure cuffs with Velcro closure and D-ring:

NOTICE! Damage to device

Never autoclave the cuffs.

Washing	Once the bladder has been removed, at up to 30 °C + detergent for delicates.
	Fold covers with Velcro closures before washing so that the two halves of the Velcro tape are in contact.
Cleaning by spraying or wiping with disinfectant	<p>We recommend the following cleaning products and disinfectants:</p> <ul style="list-style-type: none"> • mild soap • Bacillol AF • Meliseptol • Oxivir Plus spray • Terralin Protect • 70%-strength isopropyl alcohol • Fugaten® spray (Lysoform spray) <p>Following cleaning/disinfecting, wipe over with clean water to remove all residues and then leave to air-dry. Please use only soft cotton or microfiber cloths.</p>

5. SERVICING

5.1 Servicing the pressure gages

The pressure gage requires no servicing.

5.2 Calibration

seca recommends having the device recalibrated every two years. This can be performed by seca Service or by authorized maintenance services. You can contact the seca Service team at service@seca.com.


NOTE

- Please note that on the seca b10, b11 and b12 models, the head of the pressure gage can be conveniently replaced with a new one instead of having the metrological check [MTK] performed.
- The seca b10, b11 and b12 pressure gages are not designed to be calibrated.
- Replacing pressure gage head on seca b10, b11 and b12:

1.

2.

3.

4.

5.

6.

6. DISPOSAL

This device consists of high-quality raw materials which can be recycled. Dispose of the device in accordance with the local laws and regulations applicable in your country.

7. WARRANTY

There is a two-year warranty period from delivery for defects attributable to poor materials or workmanship. Defects which come under the warranty will be repaired for the customer free of charge against proof of purchase. Additional claims cannot be considered. Costs of transport to and from seca are the responsibility of the customer if the device is located somewhere other than the customer's headquarters. In the event of transport damage, warranty claims can only be made if the complete original packaging was used for transport and the device was secured and fastened therein according to its originally-packed condition. You should therefore keep all packaging parts.

The warranty will be voided if the device is opened by persons not expressly authorized by seca to do so.

We ask customers overseas to contact the seller in their respective country directly in the event of warranty claims.

8. DECLARATION OF CONFORMITY

This product complies with Council Directive 93/42/EEC pertaining to medical devices. For more information on the conformity of seca products, please go to www.seca.com/products.

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